

Small.
Smart.
Strong.

Introducing Tru-D iQ Scout

Tru-D iQ
A PDI SOLUTION



PDI

Big UVC power for small spaces.

Tru-D iQ Scout is a compact, mobile UVC-dosing device with the power to solve disinfection needs of small spaces.



Weighing less than 30 pounds, Scout is ergonomically designed for versatility and allows for UVC disinfection of areas such as bathrooms, exam rooms and ambulances.

Scout connects seamlessly to the **Tru-D iQ Controller**, which is the mobile command center for Scout. For easy transportation, Scout docks onto the Controller, and its lamps are protected in the base during transportation and manually deployed when ready for use.

Smart.

- Sensor360® technology ensures thorough room disinfection, every time
- Eight UVC sensors in a 360° array determine precise UVC cycle time
- PIR sensors recognize heat signature of human motion and/or presence in the space and prevent UVC from being emitted
- Connects to the **Tru-D SmartHUB** reporting system, an all-inclusive, web-based application
- Effectively inactivates unwanted microorganisms in the space

Effective.

Continuing **Tru-D SmartUVC's** commitment to and tradition of efficacy, Scout inactivates *C. diff*, *C. auris*, MSSA, SARS-CoV-2 and influenza A in both direct and shadowed areas.

Tru-D iQ
The evolution of intelligence



The Tru-D family of products.

Scout is part of the **Tru-D** iQ portfolio of products, which includes the **Tru-D** iQ Tower, the most powerful device available from **Tru-D** SmartUVC.



800-774-5799
info@tru-d.com



Tru-D[®]



TRU-D.COM

PDIHC.COM

©2023 Tru-D
PDI01238533